


X-Ray Protection Clothing

- Option Line -



Solutions

The **Option Line** for X-Ray Protection Clothing completes now our range of X-Ray Protection products.

The Option Line allows us to offer new youthful and varied aprons, shields and other protective accessoires accompanied by high-technology and continued progresses.

The Option Line impresses with the soft cover material, which is antimicrobial and water-resistant. The smooth surface is easy to clean and wrinkle free.

The Option Line is available in 0.25 / 0.35 and 0.50 mm Pb.

The Option Line is available in three different inner materials:

LLU – LiteLead / lead (standard, if not indicated)

MLU – MaxLite / light version, lead

LGU – LiteGreen / lead free

Special versions like “Frill“ and “PC material“ are optional to order. More details, info and prices upon request.

Option Line - Overview

Available in:

Coat Apron /	Coat Apron Junior (0.50 mm Pb)
DS Apron –	Maxx or Optixx Version
DS Apron –	Junior (0.50 mm Pb)
Vest –	Maxx or Optixx Version
Skirt	Dental Apron (0.35 mm Pb)
Patient Half Apron	Surgical D-O Apron
Thyroid Shield /	Gonad Shield (1.00 mm Pb)
Ovarian Shield (1.00 mm Pb)	Head Protection (0.25 mm Pb)
Patient Shield (add. 1.00 mm Pb)	Panoramic Apron (0.35 mm Pb)



With press buttons for storage or an easy mounting.



Fastening: Velcro fastening with extra strong velcro fastener.
(Skirt: Combo of Velcro and Clip)



Standard colour:

Royal blue and
if applicable Camel contrast ribbon

Other colours, colourful (mixed) versions, fastenings and inner material are possible and need only to be indicated in case of your order

Available colours:



Business Contact:

CAWO Solutions
P.O. Box 11 29
D - 86521 Schrobenhausen
Phone: (++49) 8252 9109-21
Internet: www.cawo.com
E-mail: cawo@agfa.com

Coat Apron

Option Line

Available in:

0.25 / 0.35 / 0.50 mm Pb

“S” - small

55 cm x 90/100/110/120 cm

“M” - medium

60 cm x 100/110/120/130 cm

“L” - large

65 cm x 100/110/120/130 cm



Fastening:

Velcro fastening with extra strong velcro fastener



Front pocket:

With press button for an easy mounting of a name ticker, (also afterwards)

Skirt and Vest

Option Line

Available in:

Front: 0.25 / 0.35 / 0.50 mm Pb
back: 0.25 mm Pb

Vest:

Maxx and Optixx –
two front closing versions

“S” – small

50 cm x 50/57/63/68 cm

“M” - medium

55 cm x 57/63/68/72 cm

“L” - large

60 cm x 57/63/68/72 cm

“XL” - Xlarge

65 cm x 57/63/68/72 cm

“XXL” - XXlarge

70 cm x 57/63/68/72 cm

Skirt:

“S” - small

42 cm x 50/53/57/62 cm

“M” - medium

44 cm x 53/57/62/68 cm

“L” - large

48 cm x 53/57/62/68 cm

“XL” - Xlarge

54 cm x 53/57/62/68 cm

“XXL” - XXlarge

60 cm x 53/57/62/68 cm

Business Contact:

CAWO Solutions

P.O. Box 11 29

D - 86521 Schrobenhausen

Phone: (++49) 8252 9109-21

Internet: www.cawo.com

E-mail: cawo@agfa.com

Optixx Vest



Maxx Vest



Fastening:

Vest: Extra strong Velcro

Skirt: Combo of Velcro & Clip

Skirt:

The short back and the bell shape
(smaller on the waist, wider below)
allows a more comfortable wearing.



 **Solutions**

DS Apron

Option Line

Front: 0.25 / 0.35 / 0.50 mm Pb
back: 0.25 mm Pb

Available in:
Maxx and Optixx –
two front closing versions

“S” - small
50 cm x 90/100/110/120 cm

“M” - medium
55 cm x 100/110/120/130 cm

“L” - large
60 cm x 100/110/120/130 cm

“XL” - Xlarge
65 cm x 100/110/120/130 cm

“XXL” - XXlarge
70 cm x 100/110/120/130 cm

Optixx
DS Apron



Maxx
DS Apron



Belt:

The belt enables a better wearing comfort. Press buttons sidewise and on the shoulder for storage belt, cap or collar.



Fastening:

Velcro fastening with extra strong velcro fastener.

Short back:

for a more comfortable wearing.

