

Hygienic Bags

These bags are typically used to protect digital imaging plates in intra-oral dental radiology applications.



Hygienic Bags are available in following sizes:

- Dental size 0 (46 x 32 mm \pm 1mm)
- Dental size 1 (51 x 34 mm \pm 1mm)
- Dental size 2 (52 x 41 mm \pm 1mm)
- Dental size 3 (65 x 37 mm \pm 1mm)
- Dental size 4 (88 x 68 mm \pm 1mm)

The minimum order quantity of each size = 5 boxes
(1 box contains 500 bags)



Business Contact:
CAWO Solutions
P.O. Box 11 29
D - 86521 Schrobenhausen
Phone: (++49) 8252 9109-21
Internet: www.cawo.com
E-mail: cawo@agfa.com



Solutions