

CAWO Gridholder

The CAWO Grid Holder can be applied with CAWO ABS and MMS-FA Cassettes (with or without window)

It is used for two reasons:

- The required grid is too thick to be mounted into the cassette
- One grid is used for more cassettes of the same size

The grid is mounted between the two plexiglas plates. The cassette is pushed into the holder and the clips in the frame keep the cassette safely in the holder



Available Sizes

- 18 x 24 cm
- 18 x 43 cm
- 24 x 30 cm
- 30 x 40 cm
- 35 x 35 cm
- 35 x 43 cm

Business Contact:
CAWO Solutions
P.O. Box 11 29
D - 86521 Schrobenhausen
Phone: (++49) 8252 9109-21
Internet: www.cawo.com
E-mail: cawo@agfa.com



Solutions