

# SAFELIGHT LED-X

LED safelight ( $\lambda=640$  nm) for use in darkrooms for processing all blue and orthochrom. medical x-ray films.



## Technical specifications:

<i>Electrical Power:</i>	AC-DC adapter: AC 110-230V / 50-60Hz DC 12V / 700mA
<i>Wattage:</i>	approx. 5 W
<i>Dimensions (H x W x D):</i>	approx. 16 x 11 x 7 cm
<i>Weight:</i>	approx. 0.4 kg
<i>Material:</i>	Chassis: ABS white Filter: PMMA
<i>Filters:</i>	red LED and PMMA-filter ( $\lambda=640$ nm)
<i>Cable length:</i>	approx. 1.5m
<i>Operating conditions:</i>	Temperature: -30°C ...+50°C Relative humidity: 30% ... 90%

**Business Contact:**  
CAWO Solutions  
P.O. Box 11 29  
D - 86521 Schrobenhausen  
Phone: (++49) 8252 9109-21  
Internet: [www.cawo.com](http://www.cawo.com)  
E-mail: [cawo@agfa.com](mailto:cawo@agfa.com)

 **Solutions**