

# Darkroom Accessories

## Processing tank set - 3 tanks -

For Film Hangers up to size  
35 x 43 (40 x 40 cm),  
with a top bar 46 cm consisting of:

- 2 tanks, each 20 liters with lid and drain cock for fixer and developer
- 1 tank, 30 ltr. with inlet and overflow, for intermediale and final washing without lid
- 1 connection set:  
2 m flexible tube Ø 25 mm,  
1 hose clamp
- 1 drip tray with drain system

material: made in pvc dark gray  
lid: made in pvc red  
space required: 50 x 60 x 62 cm  
weight: 22,00 kg

capacities:

- Developer 20 l
- fixer 20 l
- water 30 l



CAWO Version 2018-01

### Business Contact:

CAWO Solutions

P.O. Box 11 29

D - 86521 Schrobenhausen

Phone: (++49) 8252 9109-21

Internet: [www.cawo.com](http://www.cawo.com)

E-mail: [cawo@agfa.com](mailto:cawo@agfa.com)



# Solutions