

X-Ray Accessories

Cassette Holder Wall Mounted

Metal white powder coated, mounted to the wall by four screws.

The cassette holding channels slide freely along the vertical rails.

Cassettes either in vertical or horizontal position.

Suitable for any size of cassettes

dimension: 103 x 54 x 14 cm

net weight: 7.50 kg



Cassette Holder Standing

For all cassettes format with tilt-free, running liner, height adjustable with interlock button. vertical rails and cassette holder - stainless steel, cross braces and liners - made of metal with white powder coating

Dimensions: 49 x 56 x 206 cm

Net weight: 22 kg

Also available, upon request:
-mobile with 4 antistatic locking castors

- handles (on both sides)



CAWO Version 2018-01

Business Contact:

CAWO Solutions

P.O. Box 11 29

D - 86521 Schrobenhausen

Phone: (++49) 8252 9109-21

Internet: www.cawo.com

E-mail: cawo@agfa.com



Solutions