




Darkroom Accessories

Cassette Hatch

The Cassette Hatch enables you to transfer cassettes from the X-ray room into the darkroom saving the darkroom from the incidence of daylight and X-ray radiation.

Technical Data:

- Housing: powder-coated steel
- Bottom: specially coated for safety of cassettes
- Doors: 180 ° opening
lightproof interlocking



Business Contact:
CAWO Solutions
P.O. Box 11 29
D - 86521 Schrobenhausen
Phone: (++49) 8252 9109-21
Internet: www.cawo.com
E-mail: cawo@agfa.com