

Film Markers



CAWO Version 2021-04

Lead Markers are used to expose marks and legends on the X-ray film during exposure. They are transparent plastic tablets engraved with the letter or mark which is filled up with a radiopaque lead plastic. They are about 3 mm thick.

Letters and numbers: 7,5 or 10 mm high
R & L clips

L Markers: rectangular 15 mm high
R & L on chain

all Film Markers are *subject to prior sale*

Business Contact:
CAWO Solutions
P.O. Box 11 29
D - 86521 Schrobenhausen
Phone: (++49) 8252 9109-21
Internet: www.cawo.com
E-mail: cawo@agfa.com

 **Solutions**