


# Viewing Boxes

## BASic series



a comfortable support for examining  
your X-ray photographs



# Solutions

## Illuminators, BASic Series

Economical: this range, made of ABS, provides you with a comfortable support for examining your X-ray photographs.

Medical diagnosis device (class 1).

## Technical Characteristics

- White box frame: made of extruded twin-materials PVC and injected ABS. 10/10 sheet steel rear plate with epoxy coating
- Wall fixing: using fixing slots at the back of the unit.
- Light source: fluorescent tubes, “daylight” type, colour temperature >5500°Kelvin.
- Diffuser screen: Perspex, thickness 3 mm.
- Film holder: steel rollers. Hardwearing and totally reliable.
- Lighting control: bipolar switch (I) with incorporated indicator light.
- Viewing areas: extended only
- Mains cord: standard, 2 poles + earth.
- Mains electricity supply: 230V, 50 Hz (for other voltages, consult us).

CAWO Version 2020-03

### Specifications

Type	lighting sections	Viewing area	Overall measurements			Net weight	Number of tubes
			Length	Height	Depth		
horizontal BASic N 40	1	36 x 43 cm	39,5cm	51,0cm	11,0cm	5kg	1 x 32W
horizontal BASic N80	1	72 x 43 cm	75,5cm	51,0cm	11,0cm	10kg	2 x 32W
horizontal BASic N120	1	121 x 43 cm	125,5cm	51,0cm	11,0cm	13kg	2 x 36W
horizontal BASic N150	1	151 x 43 cm	154,5cm	51,0cm	11,0cm	16kg	2 x 58W

### Packaging

Type	Measurements	Volume	Weight
BASic N40	50 x 62 x 18cm	0,056m <sup>3</sup>	7kg
BASic N80	87 x 62 x 18cm	0,100m <sup>3</sup>	10kg
BASic N120	136 x 62 x 18cm	0,150m <sup>3</sup>	16kg
BASic N150	165 x 62 x 18cm	0,185m <sup>3</sup>	19kg

Business Contact:  
 CAWO Solutions  
 P.O. Box 11 29  
 D - 86521 Schrobenhausen  
 Phone: (++49) 8252 9109-21  
 Internet: [www.cawo.com](http://www.cawo.com)  
 E-mail: [cawo@agfa.com](mailto:cawo@agfa.com)

