

# X-Ray Protective Curtain

## Protective Curtains

### Lead Equivalency:

0.25, 0.35 or 0.50 mm Pb

### Covering Material:

robust and non-irritant nylon

### Colours:

see colours chart

### Dimensions:


width: 120 cm

length: 200 cm (standard)

### Protective curtain accessories:



Hanger for ceiling or for wall



Rail



Stopper



Roll



Hook



### Business Contact:

CAWO Solutions

P.O. Box 11 29

D - 86521 Schrobenhausen

Phone: (++49) 8252 9109-21

Internet: [www.cawo.com](http://www.cawo.com)

E-mail: [cawo@agfa.com](mailto:cawo@agfa.com)

 **Solutions**